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
specific for instruments

concentrated powder | use at 2%

- Bactericidal
- Fungicidal
- Mycobactericidal, including tuberculocidal
- Virucidal

1200 g | disinfectant for medical devices | www.zhermack.com

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EN - High-level disinfectant and tri-enzymatic cleaner in powder form for ordinary, surgical and rotating instruments. The product is not suitable for objects made of copper, brass, aluminium alloys, nickel alloys or damaged nickel-chromium plated instruments.

Use

For every litre of water, add 20 g (2 scoops) of Zeta 2 Enzyme. Stir to facilitate dissolution of the powder. Wait for activation of the solution: 15 minutes with water temperature of about 35°C (warm water), at least one hour with water temperature below 35°C (room temperature water). Soak the instruments for 10 minutes. Remove the instruments from the bath, rinse and dry them. If necessary, remove stubborn stains using a brush. We recommend following the recommended soak time. The prepared solution is stable for at least 8 hours.

Range of action:

Micro organisms	EU standards	Concentration	Reaction time (minutes)
Bactericidal: (S. aureus, P. aeruginosa, E. hirae)	EN 13727	0.5%	5
	EN 14561	0.5%	5
Fungicide resistance: (C. albicans)	EN 13624	0.5%	5
	EN 14562	0.5%	5
Mycobactericidal, including tuberculocidal (M. terrae, M. avium)	EN 14348	2%	10
	EN 14563	2%	10
Virucidal: (tested on Polioviruses, Adenoviruses, Parvovirus including HIV, HBV, HCV)	EN 14476	2%	1

Composition

Oxygen based (> 25%) bleach; activator (> 18%), anionic and non-ionic surfactants (<5%), corrosion inhibitors, tri-enzymatic complex (protease, lipase, amylase), additives and auxiliaries.

Warnings, storage conditions and precautions

The product is intended for use by qualified operators, with knowledge of the required safety standards. Use protective gloves and appropriate PPE. Operate in well-ventilated areas. Keep the product tightly closed in a dry and cool place. **IMPORTANT:** the undissolved powder remaining on the bottom guarantees the effectiveness of the solution for its entire duration in accordance with the instructions for use. Precautions: rinse and clean the dilution pan before introducing Zeta 2 Enzyme. Do not mix Zeta 2 Enzyme with any other disinfectants.

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