

Type 6 Self Adhesive Indicator (3.5 min at 134°C cycle)


The Valisafe Type 6 Self Adhesive cycle verification indicators offer the highest level of proof of the achievement of Sterility Assurance levels in your trays and items.

- (i) These indicators are calibrated specifically to your cycle requirements.
- (ii) Valisafe Type 6 indicators offer validated security that every tray or item has been effectively processed, thereby providing a direct link between the sterility of the items and the patient.
- (iii) Self Adhesive to provide ease of use in any document or traceability system.

The products bear the BSI Kitemark logo, this provides assurance that samples are regularly subjected to rigorous, independent testing by the British Standards Institution (BSI) to ensure they comply with stringent standards for safety, product performance, and reliability.

- Conforms to EN ISO 11140–1 Type 6
- Self Adhesive backing provides fixture to patient records
- Specifically calibrated to 3.5mins - 134°C
- Monitors critical parameters
 - Time
 - Temperature
 - Steam
- Smaller & more compact design for use in every tray or item
- Clear and abrupt colour change
- Easy to interpret - easy to record

Name:	Type 6 Self Adhesive Indicator
Classification:	EN-ISO 11140-1 Type 6
Operating temp & time:	134°C - 3.5mins
Storage conditions:	Cool, dark, dry conditions before and after use 0°C - 30°C
RH:	30 - 60%
Correct symbols:	In accordance with EN-ISO 11140-1 (Type 6)
Indicator:	Type 6
Self Adhesive:	Yes
Disp/single use:	Yes
Lot no.	Yes on each indicator & primary box
Exp date:	Yes on each indicator & primary box
Shelf life:	3 years from date of manufacture
Indicator size:	40 ± 2mm x 40 ± 2mm
Primary Box Size:	45 x 45 x 94mm
Pass/fail result:	Yes
Colour change:	Yellow to blue/ purple after exposure to dry saturated steam cycles.



Order details:

V3504065

Type 6 Self Adhesive Indicator - 3.5min at 134°C

Pack of 200